

Glass Sample Cells


Round and Rectangular

For Protein and Tissue studies, and Imaging applications in specialized probes

Round Cells

These cells can be produced with flat (most common), round, or no bottom. Lengths are optional with some of the most common listed below. Some examples are:

OD (mm)	Length (mm)	Wall (mm)
3.9	15	≤0.5
5.0	15	≤0.5
5.0	22	≤0.5
5.0	25	≤0.5
5.0	40	≤0.5
5.0	75	≤0.5
10.0	75	≤0.5
10.0	100	≤0.5
20.0	100	≤0.5
30.0	100	1.4



Glass Calibration Cells

Flat bottom.

ID (mm)	Length (mm)	Wall (mm)
1.8	20	≤0.5
2.4	20	≤0.5

This list represents some of the items that have been made for special sampling applications. Please inquire about any sampling accessories you may need for your special applications.

Prices upon request.

Rectangular Cells

These cells can be produced with flat bottoms or open-ended. Lengths are optional. Some examples are:

Inside Dimensions (mm)	Length (mm)	Wall (mm)
2.2 x 7.2	20	<0.5
4.3 x 6.3	20	<0.5