


## QUARTZ SAMPLE TUBES FOR EPR (ESR) SPECTROSCOPY



The sample tubes listed here are manufactured from commercial grade quartz material processed for use in EPR applications. For sizes **not shown**, please request.

Catalog Number	OD	ID	Length
NE-400-1.6-120	1.6mm	1mm	120mm
NE-400-2.9-160	2.9mm	2.5mm	160mm
NE-400-3.9-180	3.9mm	3.0mm	180mm
NE-400-4-250	4mm	3mm	250mm
NE-400-5-178	5mm	4mm	178mm
NE-400-5-250	5mm	4mm	250mm

## SAMPLE TUBES FOR RAMAN SPECTROSCOPY

A thin wall 5mm OD Borosilicate glass tube, 4" long.

Packed 25 tubes per package with caps packed separately.



Catalog Number	Description
NE-Raman-5	Raman Sample Tube, 5mm